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This manual is an integral part of the product with which it is enclosed and must be conserved for future consultation. Carefully read the recommendations and instructions contained in this instruction booklet before carrying out installation, operating, adjusting or maintenance operations, since they give important information regarding the SAFETY OF USE and MAINTENANCE.

Even if this device is provided with protections, the use of equipment without following the described procedures could produce a wrong working and equipment's damages.

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Pressure: before opening the lid, check that the manometer will show that the pressure rate has decreased to 0 bar.

Don't use the autoclave in lack of water Each sterilization cycle consumes approximately 1 liter of water.

At each end of the cycle, check that there are 3.5 liters of water inside, or that the heating element and probe are covered by water.

CERTIFICATO

Nr. 50 100 3290 Rev.006

SI ATTESTA CHE / THIS IS TO CERTIFY THAT

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È CONFORME AI REQUISITI DELLA NORMA
HAS BEEN FOUND TO COMPLY WITH THE REQUIREMENTS OF

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Design, production and service of scientific appliance and laboratory equipment. Trade of laboratory equipments, disposable products, glassware and medical items for scientific and technical purposes (IAF 19, 29)

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TESTING DECLARATION

Asal declare to test the autoclave 760 NE, c/o its seat in Cernusco S/N (Mi), before the sale.

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Valvola di sicurezza certificata PED tarata a / Safety valve PED certified and set to: 2 bar

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
10 – DEACTIVATION AND DISMANTLING

Compliance with WEEE Directive

Italian legislative decree 14 marzo 2014, n. 49, implementation of WEEE Directive 2012/19/UE, rules the criteria for disposal and recovery of electrical and electronic materials. All the products sold by ASAL s.r.l., are professional equipments, belong to the category 8 annex 1 Dlgs 14 marzo 2014 n. 49, so they shall be subject to the directive. The product should not be disposed of as household waste not recyclable. Use the appropriate structures of local WEEE collection for disposal of this product or call dealer or manufacturer, subscribed to the "National Register of subjects required to finance the management systems WEEE" and the collective system for the management of WEEE called "Consorzio ECOLIGHT".

The manufacturer is obliged to retire if the equipment produced after 13 August 2005; if produced before that date, the withdrawal is subject to the purchase by the customer for an equivalent product, and used the same functions (art. 24 Dlgs 14 marzo 2014 n. 49).

Disposal inappropriate may have negative consequences for the environment and human health. The recycling of materials helps conserve natural resources. The disposal of WEEE as municipal waste is punished with sanctions specifically provided for in the case of illegal disposal of this waste.

This symbol on the product or packaging indicates that the same once finished, its use must be delivered to a collection centre for the recycling of electrical and electronic equipment.

Only trained and equipped personnel have to do these operations, and they must be in accordance with Dlgs 14 marzo 2014, n. 49 and European directive 2012/19/UE for the treatment and recycle of the electric and electronic wastes.

Disconnect the power rated (see the paragraph 3 - installation). Strip these parts:

1. Electric parts
2. Plastic parts
3. Plexiglas parts
4. Iron parts
5. Syntetic parts

Follow the instructions of one's municipality.

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file: user guide
mod. 760
page: 19
edition: 03/24

To restore operation:

- ◆ turn off from the green main switch
- ◆ leave to cool until the pointer on the pressure gauge indicates “0 atm/bar” to signal the lack of pressure inside the chamber.
- ◆ open the lid and remove the basket, pour the necessary water into the sterilization chamber to completely cover the resistance
- ◆ restart by repeating the operations listed in chapter 4 - OPERATION

7 - GUARANTEE

- ◆ The good working of the described product is guaranteed from ASAL for a period of 12 months from the purchase date.
- ◆ The guarantee includes the repair or the free exchange of the instrument components that ASAL technicians admit to be defective in terms of manufacture or material (carriage and packing costs to be paid by the buyer).
- ◆ The guarantee will be invalidated by negligence, wrong installation and uses which are not observed in this instructions book, incident and tampering or transfer damage (excluding when it is made by our authorized people).
- ◆ The guarantee extension is no longer valid after repair.
- ◆ In case of restitution, the equipment has to be returned clean in order to avoid any problem for our personnel. We will refuse to draw back all the instruments that do not observe these sanitary regulations.
- ◆ Fuses are not included in the guarantee.

8 - DOCUMENT VALIDITY

ASAL reserves the right to review this document without any notice (any variations do not prejudice the observation of the directives).

9 - SPARE PARTS

See page 22-23-24.

For any spare parts needed contact the builder (see address and telephone number on the last page).

file: user guide
mod. 760
page: 18
edition: 03/24

5 - MAINTENANCE

Cleaning:

For the box powdered painted: use domestic detergents or cleaners.

For the vessel inner surface in stainless steel it is suggested use soapy water, detergents for glass or denatured alcohol, however not-abrasive products.

Do not use products with chlorine for the stainless steel surface.

The maintenance has to be carried out by experts.

To look at the good state of the autoclave and the plug.

In order to avoid any contact with hot parts, the operations have to be done in total safety conditions (disconnect the plug from the supply power cable).

To keep the autoclave in perfect working order, we recommend the following:

- ◆ the lid sealing gasket must be kept clean (in case of replacement, make sure that it is perfectly inserted in its seat).
- ◆ use deionized or decalcified water to protect against the formation of scale deposits.
- ◆ Every 5 cycles it is advisable to replace the water to avoid the onset of limescale in the short term
- ◆ To empty the vessel, use the drain cock located on the rear side with side cock and have it emptied into a suitable container.

6 - WORKING DEFECTS

In case of lack of tension, autoclave stop to work.

At the return of tension, the heating restore, check the setted time sterilization on the timer.

In case of failure to ignite:

- ◆ check the fuses
- ◆ check that the outlet is in a good state

In the event of a lack of water above the level of the resistance, inside the sterilization chamber, a protection device with acoustic alarm intervenes and the "no water" light comes on (fig. 1 pos. 3).

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7. To switch on the appliance, use the green luminous main switch. The temperature of the water inside the sterilization vessel is displayed on the digital thermometer.
8. Set the desired sterilization time using the 0-120 minute timer bearing in mind that the time necessary for heating must be added to the sterilization time, approximately 15-20 minutes to reach 121°C corresponding to 1 atm/bar.
9. When the pressure gauge indicates, after about 15-20 minutes, that the autoclave has reached pressure (1 atm/bar - 121° C), the sterilization phase begins. Verify that the time left is at least 15 minutes.
10. At this point you will notice a steam release from the counterweight of the discharge valve (calibrated at 1 atm/bar.) placed on the lid: this serves to keep the pressure constant inside the autoclave.
11. Leave it running for the set sterilization time.
12. At the end of the sterilization cycle, the timer has returned to 0, switch off the appliance using the main switch.

IMPORTANT: Before opening the sterilization vessel, allow it to cool until the pointer on the pressure gauge indicates “0 atm/bar” to signal the lack of pressure inside the vessel.

To open the lid, unscrew the 3 handwheels with joint and rotate on the axis of the fixed left handwheel.

Each sterilization cycle consumes approximately 1 liter of water.

At each end of the cycle, check that there are at least 3.5 liters of water inside, or that the heating element and probe are covered by water.

Every 5 cycles it is advisable to replace the water to avoid the onset of limescale in the short term.

To empty the vessel, use the drain cock located on the rear side with side cock by connecting the supplied silicone tube and have it emptied into a suitable container only after the pressure gauge indicates 0 atm/bar.

ATTENTION: with a pressure gauge at 0 atm/bar, hot water at around 90°C comes out of the drain cock.

3 - INSTALLATION

In order to have a functional working position, it is necessary to put the item on a stable support. To move it lift from lower part.

Don't use in dangerous atmosphere. Don't use to heat dangerous material or toxic gases that can bring about implosion or explosion, without precaution.

Before the connection of the equipment check that:

- ◆ the rated voltage value indicated on the equipment nameplate corresponds to the actual mains values.
- ◆ the power supply has an effective connection on ground.
- ◆ the electric controls are in a not-working position.
- ◆ plug in the power supply cable.

4 - INSTRUCTION FOR USE

3.5 liters of deionized water or just enough to cover the resistance are necessary for sterilization.

IT' RECOMENDED DO NOT USE DISTILLED WATER BUT DEIONIZED OR DECALCIFIED WATER.

During all operations, wear PPE suitable for the use of a pressure equipment (anti-scald gloves, goggles, safety shoes)

1. Unscrew the cover fixing handwheels. The left steering wheel is the only one that just needs to be loosened. Once the handwheels have been released, turn the cover to the left.
2. Remove the basket and pour 3.5 liters of deionized or decalcified water into the sterilization chamber, or just enough to cover the resistance.
3. Always use the internal support for the basket, on the bottom of the boiler to avoid damaging the temperature probe and avoid direct contact with the resistance.
4. Reintroduce the basket with the objects to be sterilized.
5. Close the lid, making it adhere to the gasket on the pan flange, and close it by tightening the fixing handwheels.
6. Insert the drain valve (pos. 11 fig.1 page 22) in the outlet nozzle on the lid, pressing lightly until you hear it click into place.

TABLE OF CONTENTS

1	DESCRIPTION	page 14
2	OPERATING FEATURES	page 15
3	INSTALLATION	page 16
4	INSTRUCTION FOR USE	page 16
5	MAINTENANCE	page 18
6	WORKING DEFECTS	page 18
7	GUARANTEE	page 19
8	DOCUMENT VALIDITY	page 19
9	SPARE PARTS	page 19
10	DEACTIVATION AND DISMANTLING – Compliance with WEEE Directive 2012/19/EC D.lgs 14 marzo 2014, n. 49	page 20
11	WIRING DIAGRAM	page 21
12	PARTS LIST	page 22-23-24
13	ACCESSORIES	page 25
14	EU DECLARATION OF CONFORMITY	page 26
15	TABLE OF SYMBOLS AND TESTING DECLARATION	page 27

file: user guide
mod. 760
page: 14
edition: 03/24

1 - DESCRIPTION

The autoclave mod. 760 NE is a vapour sterilizer that operate at pressure 1 atm/bar, that correspond at the sterilization temperature of 121 °C.

It's equipped with:

- ◆ Stainless steel sterilization vessel, volume 17.5 litres, in accordance with the requirements of Directive 2014/68/UE (PED). Sterilization vessel made of stainless steel X5CrNi 18-10 AISI 304, internal surface electro polished.
- ◆ Fire and flame retardant vessel insulation material.
- ◆ Lid in stainless steel with mechanical locking on the chamber by 4 handwheels.
- ◆ Lid opening on rotation pin with translational movement.
- ◆ Manual water filling. 3.5 liters of preferably deionized water are necessary for sterilization.
- ◆ Water outlet on the rear side with tap and hose connection.
- ◆ Immersion electric resistance with protective support.
- ◆ Operating pressure: 1 atm/bar
- ◆ Design pressure: 2.5 atm/bar
- ◆ Hydraulic test pressure: 5 atm/bar
- ◆ PED certified safety valve calibrated at: 2 atm/bar
- ◆ Sterilization/use temperature: 121 (± 1 °C)
- ◆ Pressure gauge for checking the operating pressure: measuring range 0-4 atm/bar. Resolution ± 0.2 atm/bar.
- ◆ Safety thermostat and acoustic alarm, in case of water shortage.
- ◆ Digital thermometer: measurement range from ambient T °C to +150°C. Accuracy ± 1 °C, Resolution ± 1 °C.
- ◆ PT 100 temperature detection probe.
- ◆ Timer 0-120 min. commands the resistance.
- ◆ Main switch and indicator lights for signaling power supply on, lack of water and electric resistance on.
- ◆ 1/8" connection (with cap) for inserting a second temperature probe for execution of certified calibrations and checks.
- ◆ External casing in steel with anti-acid powder coating.
- ◆ 4 feet in non-slip rubber.
- ◆ Exempt from ISPE SL control because it is within the foreseen capacity and pressure limits.

file: user guide
mod. 760
page: 15
edition: 03/24

2 - OPERATING FEATURES

2.1 - TECHNICAL SPECIFICATIONS

- ◆ Sterilization vessel internal dimensions \emptyset x height: \emptyset 248 x 365 mm. (useful internal height 320 mm).
- ◆ Vessel Capacity: 17.5 litres
- ◆ Overall dimensions WxDxH: 375 x 390 x 525 mm.
- ◆ Weight: 22 Kg.

2.2 - ELECTRICAL SPECIFICATION

- ◆ Power supply: 230 V - 50 Hz.
- ◆ Electrical input: 1500 W
- ◆ Overload fuses: 2 x 16 AF (5 x 20) mm.
- ◆ Connection outlet: 16 A

2.3 - CLIMATIC CONDITION FOR STORAGE, PACKING AND USE

- ◆ Altitude up to 2000 m
- ◆ Temperature 5 °C to 40 °C
- ◆ Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40 °C
- ◆ Mains supply voltage fluctuations not to exceed $\pm 10\%$ of nominal voltage
- ◆ Transient overvoltages II
- ◆ Pollution degree 3