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LyoQuest

INDEX / ÍNDICE

1. GENERAL DATA	6
2. DETAILS OF THE FREEZE-DRYER	7
2.1. Aim of this manual	7
2.2. Application	7
2.3. Requirements	7
2.4. Description	7
3. TECHNICAL CHARACTERISTICS	9
4. TRANSPORT AND PACKING	10
5. BEFORE START-UP	10
5.1. Placing on site	10
5.2. Electrical connection	10
5.3. Vacuum pump	10
5.4. Refrigerating unit	11
6. RECOMMENDATIONS	11
6.1. Ice condenser	11
6.2. Greasing	11
6.3. Long stoppage	11
7. WORKING PRINCIPLE	12
7.1. Freeze	12
7.2. Lyophilization	12
7.2.1. <i>In manifold</i>	12
7.2.2. <i>In chamber</i>	13
7.2.3. <i>End of the process</i>	14
7.3. Lyophilization modes	14
7.4. Selectable processes	14
8. OPTIONAL ACCESORIES (COMPILATION AT THE END OF THE MANUAL)	15
8.1. Chambers	15
8.2. Manifolds	16
8.3. Flasks condensations collector (Drip tray)	16
8.4. Lyoquestlogger software	17
9. INSTRUCTIONS FOR USE	18
9.1. Start up	18
9.2. Semiautomatic mode	19
9.3. Automatic mode	19
9.4. Defrosting mode	20
9.5. Menu	20
9.6. Hours counter vacuum pump	21
9.7. Stop of the freeze-dryer	21
9.8. No Plus version	22
10. RECIPES	23
11. ALARMS	23
12. MAINTENANCE PLAN	25
12.1. Vacuum unit maintenance	25
12.2. Refrigerating system maintenance	27
12.3. Surface cleaning and disinfection	27
12.4. Maintenance table	27
13. TROUBLESHOOTING	28
13.1. Vacuum unit	28
13.2. Refrigerating unit	28
13.3. Freeze-dryer kit	29

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1. GENERAL DATA

This instruction manual is valid for the laboratory freeze-dryers models:

LYOQUEST -55 (PLUS / NO PLUS)

LYOQUEST -85 (PLUS / NO PLUS)

CAUTION


- The reading and comprehension of this manual are mandatory before the start up of the freeze-dryer.
- All safety instructions that are mentioned in this manual must be respected.
- Before any maintenance operations can begin the freeze-dryer must be disconnected from the power supply.
- A use or an application not according to the recommendations expressed in the manual can generate a risk situation and means the loss of the manufacture guarantee.

The European Union issued the Directive 2012/19/EU on Waste of Electrical and Electronic Equipment (WEEE).

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive. Mark symbol is shown at the right.

Azbil Telstar sells products through distributors throughout Europe. Contact your local distributor for recycling/disposal.

The name plate & CE mark sticker is located at the back of the freeze-drier.

MANUFACTURER: AZBIL TELSTAR TECHNOLOGIES, S.L.U.

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2. DETAILS OF THE FREEZE-DRYER

2.1. AIM OF THIS MANUAL

The main aim of this manual is to act as a guideline for the handling, operation and maintenance of the equipments LYOQUEST -55 and LYOQUEST -85.

2.2. APPLICATION

This model of freeze-dryer can be used for the freeze-drying in flasks, vials, ampoules or in bulk, by means of the different configurations of the basic unit and the possibility of coupling several accessories.

WARNING

- This freeze dryer is NOT a sterilizable unit.
- This freeze dryer consists of methacrylate parts. Not following the instructions about the disinfection of the unit (chapter 12.3) or in the event of freeze drying methacrylate corrosive products can damage the parts, being excluded of the guarantee of the unit.
- Under no circumstances may the unit be used with corrosive, flammable or explosive gases.
- In the event of freeze drying products containing corrosive or explosive solvents, please contact the Customer Service, to take the specific preventive actions for the specific product.
- Certain freeze-drying cycles increase the required preventive maintenance actions on the vacuum pump. See chapter 12.1.

2.3. REQUIREMENTS

Electrical supply: 230V monophasic + PE, 50 Hz / 230V monophasic +PE, 60 Hz

Installed power: LYOQUEST -55: 2.5 kW
 LYOQUEST -85: 3.5 kW

Consumption in amperes: LYOQUEST -55: 7A (11A max)
 LYOQUEST -85: 9A (15A max)

Maximum environmental temperature to work: 25°C

2.4. DESCRIPTION

The LYOQUEST freeze-dryer is constructed as a compact desktop version. The basic unit consists of a metallic box with lateral access doors to perform maintenance tasks. The vacuum pump is supplied separately and is optional. The control panel is located at the front.

LYOQUEST components:

- A two-stage vacuum pump provided with a gas ballast or a bubbler (gas purging apparatus). When the pump is switched off a closing valve isolates the pump from the freeze-dryer, preserving the latter's vacuum while atmospheric pressure is re-established in the interior of the vacuum pump, so backlash of oil is avoided. An exhaust filter is also included.
- A refrigerating system, which cools the condenser coil to the working temperature. The condensation surface has been dimensioned ensuring the condensation of the water vapor released by the products that are being freeze dried. This prevents vapors from entering the vacuum pump. The refrigerating system is formed by:
 - Compressor:
 - LYOQUEST-55: one-stage hermetic compressor, condensed by air.
 - LYOQUEST-85: two one-stage hermetic compressors, condensed by air.
 - Condenser:
 - Material: AISI 316L stainless steel.

LYOQUEST models:

- **PLUS:** Prepared to work with all accessories and heated plates.
- **NO PLUS:** Prepared to work with all accessories, but not with heated plates.

- 1. Manifold with 8 inlets
- 2. Touch screen
- 3. Draining
- 4. Electrical connection and switch
- 5. Connection RS232 for Lyologger
- 6. Vacuum pump
- 7. Electrical connection of the vacuum pump
- 8. Chamber with heated plates


Fig. 1

3. TECHNICAL CHARACTERISTICS

VACUUM PUMP (optional)		ULVAC GLD 136C	ILMVAC P6Z-101
Dimensions (length x width x height)	mm	487.5x170x249.5	500 x 300 x 340
Weight	kg	27	25.4
Suction flow (50Hz/60Hz)	m ³ /h	8.1/9.72	5.8/6.6
Vacuum at the pump inlet	mbar	See user manual of the pump	
Motor velocity	r.p.m.	-	1500
Motor power	kW	0.4	0.4
Oil	Capacity	liters	1
	Type (standard applications)	AV-30	LABOVAC 10

ICE CONDENSER		LYOQUEST -55	LYOQUEST -85
Condensation surface	m ²	0.25	
Minimum temperature (T _{room} ≤ +21°C)	°C	-55	-85
Ice capacity	kg	8 (4 kg/24 h)	8 (6 kg/24 h)
Material		AISI – 316L	

REFRIGERATING UNIT		LYOQUEST -55	LYOQUEST -85
Number of compressors		1	2
hermetic type		One stage	
Refrigeration power (at -35°C)	W	495	800
Refrigerant (free of CFCs)		R-507	R-507/ R-23
Refrigeration of the unit		By air circulation	

EQUIPMENT		LYOQUEST -55	LYOQUEST -85
Electrical connection	V, ph, Hz	230, 1, 50 / 230, 1, 60	
Total installed power	kW	2.5	3.5
Process vacuum. Unloaded, T _{room} ≤ 21°C, condenser temperature: -55°C (LYOQUEST -55), condenser temperature: -85°C (LYOQUEST -85)	mbar	< 5 x 10 ⁻²	< 3 x 10 ⁻²
Dimensions (height x width x depth)	mm	480 x 425 x 575	
Vacuum gauge		DIGITAL	
PT100 Temperature probe		DIGITAL	
Weight	kg	65	70

4. TRANSPORT AND PACKING

The freeze-dryer is delivered together with the optional accessories and documentation in a single pack; it is packed in a suitable way to prevent it from being damaged during transportation. The freeze-dryer and the accessories are also packed in an individual package.

The pack includes a tube of silicone grease, to carry-out the maintenance of the equipment seals.

The unit has to be shipped in vertical position according to the crate indications (adhesive labels).

After having unpacked the freeze-dryer, it is necessary to verify that the pack includes all the optional accessories shown in the order.

5. BEFORE START-UP

Before proceeding with the start-up of the freeze-dryer the following instructions must be followed step by step.

5.1. PLACING ON SITE

The unit must be placed on a table that has a smooth and horizontal surface. The basic unit must be firmly settled to prevent vibrations.

WARNING

- When assembling the chamber to the Lyoquest, do not forget to assembly the methacrylate base. It is required for proper working when chamber is installed.

5.2. ELECTRICAL CONNECTION

WARNING

- Before switching the Freeze-Dryer on, check that the supply voltage is the same as the one mentioned on the specification plate.
- Installation in the place where the Freeze-Dryer is connected to must meet the safety regulations in force and be according to the local regulations of the country or customer facilities.

Connect the cable of the vacuum pump to the basic unit.

Connect the electrical cable of the basic unit to the power supply (the plug wire is supplied with a European plug).

5.3. VACUUM PUMP

ULVAC GLD-136C / ILMVAC P6Z-101 vacuum pumps (Fig. 2 & 3):

The pump is supplied with the oil loaded. Verify the oil level (see section Vacuum unit maintenance).

In order to prevent the non-condensable vapors of the product of contaminating the oil of the pump, the system has Gas-Ballast function (see section Vacuum unit maintenance).

The vacuum pump is equipped with an exhaust filter, whose function is to retain oil vapors.

The electrical connection of the vacuum pump has to be connected to the base placed at the side of the unit. Furthermore, the vacuum pump admission has to be connected to the basic unit through the supplied pipe.

Connection of the vacuum pump expulsion:

To direct the vapors outside of the pump, connect the pipe to vacuum pump exhaust, removing the installed filter.

ULVAC GLD-136C


Fig. 2

ILMVAC P6Z-101

Fig. 3

5.4. REFRIGERATING UNIT

In order to assure a suitable refrigeration of the refrigerating unit (refrigerated via forced air circulation by means of a fan), all vents must be separated by at least 20cm from other equipment or the walls.

6. RECOMMENDATIONS

6.1. ICE CONDENSER

Before each operation it is necessary to open the water drain valve of the condenser in order to assure that the ice condensed during the previous operation has been totally eliminated. If water residues remain in the condenser or the drain pipe, they will be frozen as a result of the vacuum, thus possibly compromising work conditions.

6.2. GREASING

Slightly grease the airtight seals with silicone grease. The contact surfaces should always be kept clean and greased. Avoid dirt accumulation.

6.3. LONG STOPPAGE

If the freeze-dryer is stopped for a long period of time, the oil of the vacuum pump must be changed to ensure that there is no corrosive waste resulting in damage. Disconnect the unit from the power supply. Clean all the components of the equipment (manifold, chamber...), check that the condenser has been completely emptied and grease the airtight seals of the chamber, valves and condenser.

7. WORKING PRINCIPLE

7.1. FREEZE

The product inserted into the equipment can either be previously frozen or can be frozen in the ice condenser. To freeze it in the condenser, the product is placed into the plate stack (optional accessory), and the condenser is cooled selecting the option "freezing" on the touch screen, both in semiautomatic mode and in automatic mode.

The freeze temperature is proportional to the solidification point of product. Freezing speed influences the size of the crystals and, consequently, sublimation speed. Therefore, the most suitable freeze parameters for each product must be determined. It is advisable to freeze products at a temperature below -20°C for a period of 3-4 hours.

7.2. LYOPHILIZATION


7.2.1. IN MANIFOLD

It is used to lyophilize bulk products in glass flasks with a grinding mouth ($29/32\text{mm } \varnothing$). The products contained in the flasks (filling a quarter up to maximum half the flask) have to be previously congealed in a rotating bath of alcohol, freezing the product to the walls of the flask.

When the condenser has reached the programmed temperature, the following steps must be performed:

1. Connect the vacuum unit.
2. Connect the flasks or bottles (with previously frozen product) to the valves connected to the condenser.
3. Open the valves (see Fig. 4).

Then the vacuum is done inside flasks, and the sublimation of the product starts. A pressure of less than 10mbar is reached within 10 minutes. If this does not happen, the pump will stop and the alarm of the vacuum system will appear. In this case the valve-flask connections must be reviewed in order to assure that there is no a leak.

If it is necessary to connect additional flasks when the freeze-drying process has already begun, the following steps must be performed:

1. Close all valves.
2. Connect the new flask and open its corresponding valve.
3. Wait 10-15 seconds in order for the pressure to be restored in the unit.
4. Open the valves of the other flasks.

It is possible to provide external heat by means of an electric heat source, in order to speed up the drying process. The heat source should be at the suitable distance to apply the radiation heat that the product requires. A hot air heater directed directly to the flasks may also be used.

At the beginning of the process a layer of ice will form on the outside wall of the flask; this will gradually disappear as drying out continues. When the exterior of the flask is at ambient temperature, no ice remains in the product (the process has finished). Once the process has finished, the corresponding valves of the manifold are closed and air enters in the flasks. So they can be removed.

7.2.2. IN CHAMBER

The product may be freeze-dried inside the condenser or in the chamber:

- **Inside the condenser** (see Fig. 5): It is advisable to load the product only on the lower and middle plates. The higher plates do not receive cold as do the other plates because the vacuum supply is placed at top of the condenser. The equipment allows you to program recipes (see chapter 10) to run a complete cycle without making any manual actions.
- **Outside the condenser:** This is recommended when the quantity of product allows you to load the three plates. It is advisable to freeze the product inside the condenser, because the heated plates will be cold and heating rates from negative temperatures may be reached. Once the product is frozen, remove the plates of the condenser and place them in the chamber (see Fig. 6). Also it is possible to load a product that has been frozen in other equipment. The equipment allows you to program recipes (see chapter 10) to run a complete cycle without making any manual actions once the product is frozen and the chamber has been placed.

NOTE

- The product layer should be as thin as possible (recommended height <20mm). The time of the freeze-drying process depends for an important part from this height. It is very advisable to load the freeze-dryer uniformly and that the height of product is the same in all zones.
- The time of the process also depends on:
 - The heat provided to product (external if the manifold is used or by means of an electrical resistor placed inside the plates).
 - The maximum temperature accepted by the product.
 - The final degree of humidity desired.
 - The percentage of solid material in the product.
 - The vacuum degree obtained.
 - The total amount of product to be dried each time.
 - The temperature of the condenser.


Fig. 5


Fig. 6

7.2.3. END OF THE PROCESS

Press **STOP** to stop the pump and the refrigerating circuit. Open the air inlet of the valve on the side of the chamber. The air inlet has to be opened very slowly to avoid the formation of turbulences that will spread about dust in the flasks or plates. Next, the defrosting function is activated.

At the end of the process, the pump and the refrigerating unit will stop (the compressor stop is delayed in order to extract all gases from the circuit). Next, a connection valve of the manifold or the venting valve placed at the top cover has to be opened in order to ventilate the inside of the equipment.

It is necessary to connect a container to the outlet of the outflow pipe. Once this is done the manual outflow valve placed at the left side can be opened. Start defrosting mode (see chapter 9.4) in order to heat the surface of the condenser which allows you to extract the ice block in a few minutes. Also it is possible to heat the condenser until the entire ice block has been defrosted and all the water has been drained through the outflow valve. In both cases it is advisable to remove the chamber or the manifold to accelerate the defrosting.

Before the start of a new operation, check that the condenser has been totally defrosted, that there are no water traces in the overflow pipe, and that the overflow valve is closed again.

7.3. LYOPHILIZATION MODES

There are two possible lyophilisation modes: the AUTOMATIC mode (see chapter 9.3.) and the SEMIAUTOMATIC mode (see chapter 9.2.). The automatic mode allows you to establish the processes, the parameters and the time by means of a RECIPE (see chapter 10).

7.4. SELECTABLE PROCESSES

Freezing: It starts up the refrigerating system and cools the condenser. The cooling system requires a preparation time (5 minutes approx.) in which the temperature may increase slightly. After this the temperature will quickly decrease. This step must be longer than 30 min to make sure that the system can reach the minimum temperature.

Cool+Vacuum: It starts up the refrigerating system, cools the condenser and starts up the vacuum unit. The vacuum begins when the temperature of the condenser is lower than the established parameter (normally -40°C). The control of the vacuum level will depend on the Set Point programmed. The system will need a time longer than 20 min to reach the minimum vacuum level.

Heat plates: It starts up the electrical resistors of the plates. They warm up until the programmed temperature has been reached. Also it starts up the vacuum unit. The vacuum begins when the temperature of the condenser is lower than the established parameter (normally -40°C). The control of the vacuum level will depend on the Set Point programmed.

Stop: It stops the system.

Each one of the processes may have parameters to set. In the automatic cycle it is also necessary to set the time for each step. The parameters are the following ones:

	Semiautomatic	Automatic
Freezing	-	Time (HH:mm)
Cool + vacuum	Vacuum (mBar)	Vacuum (mBar) Time (HH:mm)
Heat plates	Vacuum (mBar) Plates temperature (°C)	Vacuum (mBar) Plates temperature (°C) Time (HH:mm)
Stop	-	-

When the equipment starts up, the compressor 1 and the refrigerating fans will start up. A few seconds later, the compressor 2 will also start up (Lyoquest -85).

Depending on the temperature of the condenser probe, the electrically-operated valve of the refrigerating circuit may take a certain time to start working. For this reason it is probable that the temperature increases some degrees initially.

The pump starts to work when the Set Point has been reached, and once it works there is a small delay in opening the vacuum valve.

8. OPTIONAL ACCESORIES (COMPILATION AT THE END OF THE MANUAL)

8.1. CHAMBERS

Material: Methacrylate (PMMA)

The vial stoppering devices under vacuum are manual and vials provided with closure are required.

Fig.	Description	Diameter	Distance between plates
7a	Cylindrical chamber equipped with 3 standard plates	220mm	70 mm standard
7b	Cylindrical chamber equipped with 3 heated plates	220mm	70 mm standard / 140 mm large spacer (LS)
8	Cylindrical chamber equipped with 3 heated plates and 8 inlets manifold	220mm	
9	Cylindrical chamber equipped with 3 heated plates, vial stoppering device and 8 inlets manifold	220mm	
10	Cylindrical chamber equipped with 3 heated plates and vial stoppering device	220mm	
11a	Big size cylindrical chamber equipped with 3 heated plates	350mm	71 mm standard / 142 mm large spacer (LS)
11b	Big size cylindrical chamber equipped with 3 heated plates and vial stoppering device	350mm	
12a	Big size cylindrical chamber equipped with 4 heated plates	350mm	71 mm standard
12b	Big size cylindrical chamber equipped with 4 heated plates and vial stoppering device	350mm	


Fig. 7a


Fig. 7b

Fig. 8

Fig. 9

Fig. 10


Fig. 11a


Fig. 11b


Fig. 12a


Fig. 12b

8.2. MANIFOLDS

Manifolds with 8-12-16 inlets

Horizontal manifold. Each one of the manifold inlets is provided with a 3 way valve which allows isolating and restoring the pressure in each one of the flasks independently. The valves allow the junction to grinding flasks by means of the conical adapter that they incorporate (external diameter: 29-32 mm; internal diameter 18 mm). By removing this adapter, tubes or flasks can be inserted directly in the valve.

Circular manifold with 8 inlets (fig 13a) and straight one with 12 inlets (fig 13b).

By joining together two 8 inlets manifolds it is possible to get 16 inlets (fig 13c).

Fig. 13a


Fig. 13b


Fig. 13c

Manifolds with 40-80 inlets

Manifold with 40 inlets for 13mm diameter pipes. By means of a distributor it is possible to get 80 inlets (fig 14b).


Fig. 14a


Fig. 14b

8.3. FLASKS CONDENSATIONS COLLECTOR (DRIP TRAY)

Cover to collect the condensations generated on the outer part of the flasks of the manifold (Fig. 15).

Fig. 15

8.4. LYOQUESTLOGGER SOFTWARE

Software for monitoring on PC (Fig. 16) that provides data for the vacuum and the temperatures of the condenser, plates and product.

The equipment is connected by a RS-232 serial port on PC and it sends information in real time.

A data packet is sent every second in order to know the freeze dryer state.

The serial port is automatically configured by the software.

The dataflow is unidirectional from the freeze dryer to the PC. It means that it is only possible to monitor, and not to control, the equipment connected to the PC.

Minimum computer requirements:

- Intel Pentium II o ADM Athlon Processor.
- 512 MB of RAM.
- Operating system: Windows XP, Windows 2000 o Windows 7.
- Fast IDE Hard Drive with more capacity than 5GB.
- SVGA Graphics card (32 bits colour at 1024 x 768 Pixels).


Fig. 16

9. INSTRUCTIONS FOR USE

9.1. START UP

NOTE

Before beginning a new operation:

- Open the drain valve in order to check that all water from the last operation has been totally eliminated. Close the valve again.
- Check that the condenser is covered by its transparent cover (manifold or chamber).
- Check that all valves of the chamber are closed correctly.

To start up the equipment, press the switch placed above the electrical connection (see Fig. 17), moving it to the position I.

The equipment will start up and the introduction screen will appear (Fig. 18).

Press on the middle of the screen to access to main screen (Fig. 19).

The main screen contains:

1. The condenser temperature value and the vacuum value. Press on this box to access the temperature probes screen (Fig. 20).
2. An indicator of working mode.
3. Four buttons, from left to right:
 - Access the automatic mode screen
 - Access the semiautomatic mode screen
 - Start defrosting mode
 - Access the menu screen

The first three buttons are used to start a working mode, for this reason they are disabled while the other mode is activated.

Fig. 17


Fig. 18


Fig. 19


Fig. 20

9.2. SEMIAUTOMATIC MODE

The second button of the main screen leads to semiautomatic mode (Fig. 21).

On the semiautomatic mode screen there is a selector to select the process to realize.

- Stop
- Freezing
- Cool+Vacuum
- Heat plates

(see the selectable processes description in chapter 7.4)

Once activated, the system will reach the preset temperature and pressure as soon as possible. Since the semi-automatic cycle does not allow to program time, the functions will not be time-dependent and the parameters will remain constant until the equipment is switched off.

The temperature set point button is used to establish the temperature that the plates will reach in the "heat plates" process.

The vacuum set point button is used to establish the vacuum value in the processes "heat plates" and "cool+vacuum".

The lower right button is used to return to the main screen.


Fig. 21

9.3. AUTOMATIC MODE

The automatic mode requires that a recipe has been added previously to the PLC (see chapter 10).

The first button of the main screen leads to automatic mode (Fig. 22).

On the automatic mode screen there are indicators: set points, step number of automatic cycle, process type and process time. On the lower part of the screen there are four buttons, from left to right:

- Start automatic cycle
- Stop automatic cycle
- Access the recipe in PLC screen
- Return to main screen

Pressing the third button, the recipe that will be followed appears (Fig. 23). The parameters of this recipe may be modified even when the cycle has begun. The process to follow in each step and the parameters of the step that is ongoing may not be modified. Use the button (arrow) located in the upper right or left corner of the screen to move from the screen Fig. 23 (steps 1 to 5 of the recipe) to the screen Fig. 23a (steps 6 to 10 of the recipe) and vice versa.

To edit step by step, the edit button has to be pressed.

The lower right button is used to return to the previous screen.


Fig. 22


Fig. 23


Fig. 23a

On the step by step edition screen (Fig. 24) the parameters of each step may be edited.

The lower left buttons are used to access the previous and subsequent step. The lower right button is used to return to the recipe in PLC screen.

Depending on the selected process in each step, the system allows to set a few parameters.

	Time (hh:mm)	Vacuum (mBar)	Plates temp. (°C)
Freezing			
Cool + Vacuum			
Heat plates			
Stop			

Required


Fig. 24

To set each one of the parameters, the corresponding box must be pressed. A numeric keypad will appear on the screen (Fig. 25).

The fixed parameters are final values of the step, that is, the values that the system tries to reach at the end of the step. The system follows a few rates taking as initial values the readings of the probes when the step starts.

In the first vacuum step, the system doesn't follow rates. The equipment will try to reach the vacuum value during the entire step.

9.4. DEFROSTING MODE

The third button of the main screen is used to activate and deactivate the defrosting mode. After realizing a freeze-drying cycle, it is possible to run a defrosting cycle in order to defrost the ice accumulated in the condenser and to extract the water by the drain (see Fig. 26). To perform the defrosting, heated gas flows by the condenser coil to increase its temperature.

The end of the defrosting mode takes place automatically, when a predetermined time has passed.

9.5. MENU

The fourth button of the main screen is used to access the menu (Fig. 27). The menu is used to program a recipe, to visualize an alarm (if an alarm is activated) and to access the adjustable parameters.

On the menu screen there are two buttons in the middle. The first of them leads to the recipes edition screen (see chapter 10). The second button leads to the alarm screen (see chapter 11). There are also four buttons on the lower area, from left to right:

- Access the settings screen
- Access the synoptic screen
- Access to vacuum pump operating hours counter
- Return to the main screen

On the settings screen (Fig. 28) the safety temperature may be modified. This screen contains boxes to select the use of the manifold or the plates. The three lower buttons, from left to right, are used to:

- Access the language selection screen
- Access the screen of screen adjustments
- Return to the main screen


Fig. 25


Fig. 27

Drain valve (opened)


Fig. 26

The language selection screen allows you to choose the language of the program (English, Spanish or French) by selecting the corresponding option.

On the screen of screen adjustments, the terminal brightness and contrast may be adjusted.


Fig. 28


Fig. 29

On the synoptic screen (Fig. 29) it is possible to visualize the main components of the equipment that are working at the given moment (they have the word ON next to the component). At the left top the evolution of the pressure may be displayed.

9.6. HOURS COUNTER VACUUM PUMP

In order to check the working hours of the vacuum compressor, press the third icon located on the inferior part of the menu screen:


Fig. 30

9.7. STOP OF THE FREEZE-DRYER




When an automatic, semiautomatic or defrosting cycle is stopped, the compressor may remain working for a few seconds in order to remove the gases from the circuit.

In case of an accidental stop, disconnection or power cut, the freeze-dryer will stop. When the power has been restored, the process continues as it was before the cut and the equipment will recover its previous state.

9.8. NO PLUS VERSION

The Lyoquest NO PLUS model is not prepared to work with heated plates. Due to these differences the NO PLUS program has some difference from the PLUS model. These differences affect basically to some screens.

These screens are shown below:

PLUS	NO PLUS
<p>SEMI-AUTOMATIC</p> 	<p>SEMI-AUTOMATIC</p> 
<p>AUTOMATIC</p> <p>Vacuum SP 0.100 mBar Shelves Temp SP: 35.0 °C Step number: 3 Process: Stop Time: 2 : 15 : 35</p> 	<p>AUTOMATIC</p> <p>Vacuum SP: 0.000 mBar Step number: 0 Process: Stop Time: 0 : 0 : 0</p> 
<p>SETTINGS</p> <p>Temperature to Allow Vacuum: -##.## °C</p> <p>Shelves <input checked="" type="radio"/> Manifold <input type="radio"/></p> <p>Lang. Screen Settings</p> 	<p>SETTINGS</p> <p>Temperature to Allow Vacuum: 0.0 °C</p> <p>Lang. Screen Settings</p> 
<p>AUTOMATIC STEP 1</p> <p>Process: HeatShelves</p> <p>Time (hh:mm): 00 : 00</p> <p>Vacuum: 0.000 mBar</p> <p>Shelves Temp: +00.0 °C</p> <p>Step 2</p> 	<p>AUTOMATIC STEP 1</p> <p>Process: Stop</p> <p>Time (hh:mm): 0 : 0</p> <p>Vacuum: 0.000 mBar</p> <p>Shelves Temp: 0.0 °C</p> <p>Step 2</p> 

10. RECIPES

This freeze-dryer has 4 recipes with 10 steps everyone:


Fig. 30a

11. ALARMS

The equipment has a set of alarms that inform the user in case of malfunctioning. If there is more than one active alarm, they will appear in a list (Fig. 31).

The alarms will be kept active until the problem has been solved, and the RESET button has been pressed. Pressing the RESET button only will have effect when the reason that has provoked the alarm has been resolved. Subsequently the alarm will be eliminated from the list.

The lower right button is used to minimize the screen, but it does not eliminate the alarm. The screen will appear again after a few seconds.


Fig. 31

High pressure compressor 1 alarm

ACTION: Refrigerating compressor 1 is stopped.

DESCRIPTION: High pressure-controller of the compressor 1 is indicating that there is an excessive pressure.

POSSIBLE CAUSES: Equipment is not receiving correct refrigeration.

Fans do not work.

Placement of the equipment does not allow a correct inlet of the air flow by the side.

High pressure compressor 2 alarm

ACTION: Refrigerating compressor 2 is stopped.

DESCRIPTION: High pressure-controller of the compressor 2 is indicating that there is an excessive pressure.

POSSIBLE CAUSES: Equipment is not receiving correct refrigeration.

The compressor 1 does not work correctly.

Fans do not work.

Placement of the equipment does not allow a correct inlet of the air flow by the side.

Condenser temperature probe alarm

ACTION: Perform vacuum is not possible, because the refrigeration of the condenser is not guaranteed.

DESCRIPTION: Probe does not give signal on expected range (4-20mA).

POSSIBLE CAUSES: Probe is damaged.

Probe signal converter is damaged.

Electrical connection among the converter and the PLC is not correct.

Vacuum probe alarm

DESCRIPTION: Probe does not give signal on expected range (1.5-8.5 VDC).

POSSIBLE CAUSES: Probe is damaged.

Electrical connection among the vacuum probe and the PLC is faulty.

Compressor (1 or 2) circuit-breaker alarm

ACTION: Compressor has stopped or cannot start up.

DESCRIPTION: Compressor circuit-breaker has blown.

POSSIBLE CAUSES: Consumption of compressor is excessive, because its functioning is anomalous.

Circuit-breaker is faulty.

Electrical connection among the circuit-breaker and the PLC is faulty.

Fans circuit-breaker alarm

ACTION: Fans have stopped or cannot start up.

DESCRIPTION: Fans circuit-breaker has blown.

POSSIBLE CAUSES: Consumption of the fans is excessive, because its functioning is anomalous.

Circuit-breaker is faulty.

Electrical connection among the circuit-breaker and the PLC is faulty.

Vacuum pump circuit-breaker alarm

ACTION: Vacuum pump has stopped or cannot start up.

DESCRIPTION: Vacuum pump circuit-breaker has blown.

POSSIBLE CAUSES: Consumption of the vacuum pump is excessive, because its functioning is anomalous.

Circuit-breaker is faulty.

Electrical connection among the circuit-breaker and the PLC is faulty.

Vacuum failure alarm

ACTION: Pump and vacuum valve have stopped.

DESCRIPTION: Correct vacuum has not been reached in 10 minutes.

POSSIBLE CAUSES: Some valve is opened (manifold, drain valve, venting valve...).

Manifold is incorrectly placed and its joint with the condenser is not tight.

Pump does not start up (verify the connection).

Vacuum valve has not opened (check the fuse).

Heaters circuit-breaker alarm (in case optional chamber is installed)

ACTION: Heaters have stopped or cannot start up.

DESCRIPTION: Heaters circuit-breaker has blown. It is placed beneath the touch screen.

POSSIBLE CAUSES: Consumption of heaters is excessive. Probably, this is due to a short circuit.

Circuit-breaker is faulty.

Connector of heaters has been disconnected for a long time.

Plates temperature probe alarm (in case optional chamber is installed)

ACTION: The use of the plates heaters is not possible.

DESCRIPTION: Probe does not give signal on expected range (4-20mA).

POSSIBLE CAUSES: Chamber probe connector has been disconnected for a long time.

Probe is damaged or probe signal converter is damaged.

Electrical connection among the converter and the PLC is not correct.

12. MAINTENANCE PLAN

This manual is a guide to the serviceable components of the freeze-dryer that require maintenance. It must be an active document that can be modified according to user requirements.

The maintenance plan consists of inspections and operations that must be performed from every day through to a period of once a year. The information will have to be compiled and detailed in order to be useful in the following maintenance operations.

CAUTION

- Before the maintenance operations begin, the freeze-dryer must be disconnected from the power supply.
- In general, a periodic inspection of the correct operation of all the unit components should be performed once a year.

12.1. VACUUM UNIT MAINTENANCE

The revision and maintenance operations that must be performed in the vacuum unit are as follows:

OIL CONDITION CONTROL

The oil condition can be observed via the sight glass. To carry out this operation, when the pump is at working temperature, switch off the unit.

In normal conditions the oil is clear and transparent. Oil draining and changing is recommended when it has a cloudy appearance (contaminated by water).

The vacuum unit has a Gas-Ballast function in order to prevent the non-condensable vapors of the product of contaminating the oil. The Ulvac 136 GLD-136C vacuum pump has external Gas-Ballast (see Fig. 32), and it is recommended to leave the Gas-Ballast opened ($\frac{1}{4}$ setting) during the first step of the freeze dryer cycle. The Ilmvac P6Z-101 vacuum pump has a bubbler (gas purging apparatus) instead of a gas ballast system.

OIL LEVEL CONTROL

The correct oil level will cover half of the sight glass. If the oil level is below this level, the necessary oil quantity must be added through the filling opening in order to reach the correct level. Before filling, the pump must be stopped.

VISUAL CONTROL OF OIL LOSS

Verify the absence of oil around the vacuum pump. The presence of oil around the pump indicates a possible leak in the joints, incorrect tightness of the caps, the oil level sight glass, or the saturation of the filter. The unit must be cleaned in order to be able to locate the origin of the leak.

OIL CHANGE

Oil changes can only be performed when the unit has run and is at working temperature. This change must be performed after the pump has worked for a minimum of 20 minutes, this to obtain more oil fluidity.

The steps to change the vacuum pump oil are the following ones:

1. Unscrew the drain cap to drain off old oil.
2. To completely drain the unit, it may be necessary to run the pump (3-5 secs) and stop.
3. Screw the drain cap.
4. Remove the filling cap.
5. Add the new oil through the filling hole until halfway the sight glass.
6. Screw the filling cap.

It is advisable to change the oil before any long periods of non-usage; it is harmful for the pump to contain dirty oil during long periods of non-usage.

PUMP BODY CLEANLINESS

The steps to clean the body of the pump are as follows:

1. Empty the used oil from the pump.
2. Refill the oil in the pump up to $\frac{1}{4}$ of the capacity.
3. Run the pump for approximately 30 seconds.
4. Empty the pump.
5. Repeat the steps 1, 2 and 3.
6. Fill the pump with oil, (up to halfway the sight glass).

OILS

Recommended Azbil Telstar Technologies oils are:

- **AV-30**, mineral oil for general applications in industry and laboratory.
- **CHEMOIL-500**, oil to be used at high temperatures.
- **WETOIL-8r**, oil to be used in aspiration of large quantity of water steam.

NOTE:

Observe the current environmental regulations about the waste and disposal of used oil.

ULVAC GLD-136C


Fig. 32

ILMVAC P6Z-101


Fig. 33

12.2. REFRIGERATING SYSTEM MAINTENANCE

The refrigerating system requires minimal maintenance. If due to accidents or breakages oil or refrigerant needs to be added, it will be necessary to get in contact with Azbil Telstar Technologies Customer Service.

To assure a proper refrigeration of the system, all the grids must be at a minimum distance of 20cm away from any other equipment and from the wall.

It is recommended to perform a periodic cleaning of the dust accumulated in the condenser batteries of the refrigerating system. To proceed, the lateral doors of the freeze dryer must be opened to clean the battery by means of a small paintbrush, and if it is possible by means of compressed air. This should be done particularly before summer.

When the switch is off, the compressor continues working for the time set by the collecting timer (2 min.) in order to inhale the gas of the evaporator. If the time has passed and the compressor has not stopped, the equipment will have to be checked, for it's not advisable that the compressor works without gas circulation.

Loads of refrigerant:

	LYOQUEST-55	LYOQUEST-85
R-23	-----	95 g
R-507	420 g	390 g

12.3. SURFACE CLEANING AND DISINFECTION

For both, unit and accessories, use sterile or non-use disables cloths, which do not shed particles or fibres.

To disinfect use cloths dampened with a disinfecting solution which do not damage or affect the paintwork, stainless steel or chamber glass. A previous clean to disinfecting can also be done with water.

To access to corners or difficult access zones, another type of tools will have to be used.

12.4. MAINTENANCE TABLE

The following table contains the advised values for a normal service of the freeze dryer. A number of very severe conditions of work can reduce the intervals expressed in the table, especially in the oil change operation of the vacuum pump.

OPERATION	PERIODICITY
VACUUM UNIT	
- Oil condition control	Before each operation
- Change of oil filter	1 year
- Oil level control	Before each operation
- Visual inspection of oil loss	6 months
- First oil change	After the first 150 hours of service
- Oil change	500 hours of service (as maximum)
- Pump body cleaning	1 year
- Check the performance of the vacuum unit, verifying the normal times to reach the vacuum throughout the unit, with the condenser cold P _{atm} to 1mbar < 5min	6 months
REFRIGERATING UNIT	
- Condenser panel cleanliness	6 months
- Check the performance of the refrigerating system, verifying the normal times to cool the condenser. LQ-85: +20 to -40 <5 min LQ-55: +20 to -30 <10 min	6 months
FREEZE DRYER	
- Joint cleaning and greasing	6 months

13. TROUBLESHOOTING

13.1. VACUUM UNIT

Anomaly: The vacuum pump doesn't start.

CAUSE	SOLUTION
The pump motor is not receiving voltage.	Check the connection of the on/off switch located upon the motor of the pump.
	Check the electrical connection terminals of the pump.
	Check the corresponding circuit-breaker.
	Check that the environmental temperature is not too high (<35°C), this might have caused the circuit breakage.
The pump's starter protector is blocked, because the condenser doesn't reach the correct temperature. Minimum temperature to switch on the pump: LQ-85 <-40°C / LQ-55 <-30°C	Check the temperature of the condenser.
Pump wheel is seized up. The motor axle doesn't turn.	Contact Azbil Telstar Technologies Customer Service.

Anomaly: Vacuum cannot be obtained.

CAUSE	SOLUTION
The motor turns but doesn't create vacuum.	Contact Azbil Telstar Technologies Customer Service.

Anomaly: The pump doesn't reach the vacuum pressure limit.

CAUSE	SOLUTION
Lack of oil.	Check the level of oil of the pump, and add oil up to the suitable quantity (please refer to paragraph 12.1.4.).
The oil is contaminated.	Empty the oil from the pump and fill it with new oil.

Anomaly: Abnormal heating of the pump.

CAUSE	SOLUTION
Lack of ventilation or environmental temperature too high (>35°C).	Place the pump in an area with better ventilation.

Anomaly: Mist leak from the emission channel.

CAUSE	SOLUTION
Gas-Ballast is opened too much.	Close the Gas-Ballast.
There is a leak.	Locate the leak.

13.2. REFRIGERATING UNIT

Anomaly: The refrigerating unit doesn't start up.

CAUSE	SOLUTION
The compressors don't receive voltage.	Check electrical connection terminals of the compressors.
	Check the circuit breaker corresponding to each compressor.

Anomaly: The condenser doesn't reach the minimum temperature.

CAUSE	SOLUTION
Refrigerant leak.	Contact Azbil Telstar Technologies Customer Service.
The refrigerating fans of the refrigerating unit don't work.	Check the circuit breaker corresponding to each fan.

Anomaly: Intermittent stoppage of the compressor.

CAUSE	SOLUTION
Insufficient refrigeration of the condenser component.	Place the freeze-dryer in an area with better ventilation.
	Check that the temperature in the room is lower than 25°C.
	Clean the area formed by the condenser component with compressed air and a paintbrush.
	Check that there is a suitable circulation of the forced air in the fan.

Anomaly: The compressor stops after alternative stoppages.

CAUSE	SOLUTION
The thermal protector of the motor of the compressor is disconnected.	The same solutions of the previous anomaly.
	Wait about 30 minutes to cool down the motor and unblock the protector, allowing again its operation.

13.3. FREEZE-DRYER KIT

Anomaly: Working vacuum cannot be reached.

CAUSE	SOLUTION
There is a leak in the equipment.	Check that the unit is watertight.
	Lubricate the watertight joints with silicone.
	Check that all the manifold keys are closed.
	Check that the condenser discharge valve is closed.
	Check there is no defrosting water of the condenser remaining from the last operation in the discharge tube.

Anomaly: The vacuum gauge indicates "gauge short-circuit" on the display (-100).

CAUSE	SOLUTION
Vacuum gauge short-circuit	Replace the gauge.

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







PLUS	NO PLUS
 <p>SEMI-AUTOMATIC</p> <p>Freezing Cool+ Vacuum</p> <p>Stop Heat Shelves</p> <p>Shelv SP 30.0 °C Vacuum SP 0.500 mBar</p>	 <p>SEMI-AUTOMATIC</p> <p>Freezing</p> <p>Stop Cool+ Vacuum</p> <p>Vacuum SP 0.000 mBar</p>
 <p>AUTOMATIC</p> <p>Vacuum SP 0.100 mBar</p> <p>Shelves Temp SP: 35.0 °C</p> <p>Step number: 3</p> <p>Process: Stop</p> <p>Time: 2 : 15 : 35</p>	 <p>AUTOMATIC</p> <p>Vacuum SP: 0.000 mBar</p> <p>Step number: 0</p> <p>Process: Stop</p> <p>Time: 0 : 0 : 0</p>
 <p>SETTINGS</p> <p>Temperature to Allow Vacuum: -##.# °C</p> <p>Shelves <input checked="" type="radio"/></p> <p>Manifold <input type="radio"/></p> <p>Lang. Screen Settings</p>	 <p>SETTINGS</p> <p>Temperature to Allow Vacuum: 0.0 °C</p> <p>Lang. Screen Settings</p>
 <p>AUTOMATIC STEP 1</p> <p>Process: HeatShelves</p> <p>Time (hh:mm): 00 : 00</p> <p>Vacuum: 0.000 mBar</p> <p>Shelves Temp: +00.0 °C</p> <p>Step 2</p>	 <p>AUTOMATIC STEP 1</p> <p>Process: Stop</p> <p>Time (hh:mm): 0 : 0</p> <p>Vacuum: 0.000 mBar</p> <p>Shelves Temp: 0.0 °C</p> <p>Step 2</p>

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Only the use of original components and accessories guarantees a correct working of the freeze-dryer.

DESCRIPCIÓN / DESCRIPTION	CÓDIGO / CODE
GRUPO DE VACÍO / VACUUM UNIT	
Bomba de vacío ULVAC GLD-136C (incluye filtro de expulsión) <i>Vacuum Pump ULVAC GLD-136C (includes exhaust filter)</i>	1517468
Filtro de expulsión bomba ULVAC GLD-136C <i>vacuum pump exhaust filter ULVAC GLD-136C</i>	622106
Major kit ULVAC GLD136A/C 100005736 <i>Service kit ULVAC GLD136A/C 100005736</i>	627385
Bomba de vacío ILMVAC P6Z-101 <i>Vacuum Pump ILMVAC P6Z-101</i>	72324
Filtro de expulsión 800160 for P6Z-101 <i>Cartridge for oil mist 800160 for P6Z-101</i>	630249
Major kit ILMVAC P6Z-101 302076-01 <i>Service kit ILMVAC P6Z-101 302076-01</i>	630248
Transductor de vacío THYRACONT <i>Vacuum transducer THYRACONT</i>	57403
Grasa de silicona <i>Silicone grease</i>	4716
Aceite AV-30 para bomba de vacío (1 litro) <i>Oil AV-30 for Vacuum Pump (1 litre)</i>	15263
SISTEMA FRIGORÍFICO / REFRIGERATING SYSTEM	
Compresor R-507 SC15CL Danfoss (-55°C & -85°C, 50 Hz) <i>Compressor R-507 SC15CL Danfoss (-55°C & -85°C, 50 Hz)</i>	58411
Compresor R-507 SC15CLX Danfoss (-55°C & -85°C, 60 Hz) <i>Compressor R-507 SC15CLX Danfoss (-55°C & -85°C, 60 Hz)</i>	62773
Compresor R-23 SC10CL Danfoss (-85°C, 50 Hz) <i>Compressor R-23 SC10CL Danfoss (-85°C, 50 Hz)</i>	58410
Compresor R-23 SC10CLX Danfoss (-85°C, 60 Hz) <i>Compressor R-23 SC10CLX Danfoss (-85°C, 60 Hz)</i>	62771
GENERAL	
Junta de estanqueidad del condensador <i>Condenser tightness gasket</i>	57481

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Fig. A1

This manual holds 3 exemplary recipes to facilitate the use of the equipment and provide the user with an initial starting point to start lyophilization.

Attention: *These recipes are merely intended as a tutorial and therefore they cannot be applied to any other product. The recipes were developed for use with the intended products only (mushrooms, skimmed milk, and a 5% mannitol solution). Evidently every other product will have its own optimal recipe which will depend on its intrinsic nature. By consequence these recipes can only be used in order to get a general idea of how to take a start with the lyophilization process. Further experimentation based on these recipes can nevertheless lead to the optimal recipe for each specific product.*

ANNEX A.1: RECIPE FOR MUSHROOMS LYOPHILIZATION

Step Nº	Process	Parameters		
		Time (hh:mm)	Pressure (mBar)	Temperature(°C)
1	Cool + vacuum	00:30	0,2	-
2	Cool + vacuum	01:00	0,2	-
3	Heat plates	01:00	0,2	30
4	Heat plates	99:00	0,2	30

Fig. A1: lyophilized mushrooms. *As can be seen the lyophilized mushrooms maintain the original volume and appearance. Noticeably this product has practically lost all its humidity, which results in a strongly reduced density compared to the original product. Drying products by lyophilization involves two major advantages: on the one hand it is possible to conserve the product during a much longer period of time, while on the other hand it entails a considerable reduction in transportation costs (due to its reduced weight). In order to maximize the conservation of the mushrooms it is recommended that once lyophilized, the mushrooms are inserted in a bag that previously has been vented with nitrogen. Once the mushrooms have been inserted, the bag needs to be closed immediately to avoid contamination by atmospheric humidity.*

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Fig. A2

Fig. A3

ANNEX A.2: RECIPE FOR THE SKIMMED MILK LYOPHILIZATION

Step Nº	Process	Parameters		
		Time (hh:mm)	Pressure (mBar)	Temperature (°C)
1	Cool + vacuum	00:30	0,2	-
2	Heat plates	01:00	0,2	20
3	Heat plates	24:00	0,2	20
4	Heat plates	00:15	0,001	35
5	Heat plates	06:00	0,001	35

Fig. A2: lyophilized skimmed milk. The lyophilized skimmed milk has the appearance of a well-formed, solid, consistent tablet with minor shrinkage. In order to obtain optimal conservation it is recommended to lyophilize the milk in vials, which are to be closed under vacuum atmosphere once the lyophilization process has come to an end. In this way contamination of the product by atmospheric humidity is avoided.

ANNEX A.3: RECIPE FOR A 5% MANNITOL SOLUTION LYOPHILIZATION

Step Nº	Process	Parameters		
		Time (hh:mm)	Pressure (mBar)	Temperature (°C)
1	Cool + vacuum	00:30	0,4	-
2	Heat plates	01:00	0,4	40
3	Heat plates	24:00	0,4	40
4	Heat plates	01:00	0,001	50
5	Heat plates	04:00	0,001	50

Fig. A3: lyophilized 5% mannitol solution. The lyophilized mannitol has the appearance of a well-formed, solid, consistent tablet without shrinkage. As with the skimmed milk, in order to maximize conservation it is recommended to lyophilize the product in vials that are subsequently closed under vacuum atmosphere once the lyophilization process has terminated.

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Note down in the attached tables the information about the maintenance operations performed in the vacuum pump and in the refrigerating compressor. It is necessary to annotate the number of hours operated before resetting the time meters.

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Symbol	Description	Refer.Numb.	Reference
-QP1	BULGIN Power inlet module	BZV49/Z0000/71	601598
-QF1	TELEMECANIQUE Motor circuit breaker	GB2-CD21	12016
-CN11	BULGIN IEC Connector	PX0675/63	62713
-Q1	TELEMECANIQUE Motor circuit breaker	GB2-CD14	12014
-Q2	TELEMECANIQUE Motor circuit breaker	GB2-CD10	12013
-Q3	TELEMECANIQUE Motor circuit breaker	GB2-CD09	57492
-Q4	TELEMECANIQUE Motor circuit breaker	GB2-CD09	26632
-Q5	TELEMECANIQUE Motor circuit breaker	GB2-CD09	26632
-TSC1	BEIJER Graphic terminal 3,5"	H-T40M-PA	74183
-PS1	MEAN WELL Power supply	MDR-60-24	62470
-AM1	MITSUBISHI RS232 Adpator module	FX3U-232ADP	59073
-CPU1	MITSUBISHI PLC	FX3UC-16MT/DSS	58606
-AI1	MITSUBISHI 4-Channel analog input module	FX3UC-4AD	58607
-R1	OMRON G2RS SPDT Relay with LED indicator	G2R-1-SN	12823
-R2	OMRON G2RS SPDT Relay with LED indicator	G2R-1-SN	12823
-R3	OMRON G2RS SPDT Relay with LED indicator	G2R-1-SN	12823
-R4	OMRON G2RS SPDT Relay with LED indicator	G2R-1-SN	12823
-R5	OMRON G2RS SPDT Relay with LED indicator	G2R-1-SN	12823
-R6	OMRON G2RS SPDT Relay with LED indicator	G2R-1-SN	12823
-R7	OMRON G2RS SPDT Relay with LED indicator	G2R-1-SN	12823
-SSR1	OMRON Solid state relay	G3NA-205B	65009
-CV1	CONVERSOR PT100 TX-ISOPACK+CAJA 65x50x45	7885882	633451
-SD1	WIKA PT-100 Sensor	WE16105	606927
-VT1	THYRACONT Vacuum transducer	VSP62MV	57403
-CV1	CONVERTIDOR PT100 -100/100°C TxIsoRail	7577317	633365
-SD2	WIKA PT-100 Sensor	WE16105	606926
-CV1	CONVERTIDOR PT100 -100/100°C TxIsoRail	7577317	633365
-SD3	WIKA PT-100 Sensor	WE16105	606926
-F1	Fuse (2 A)		2429
-F2	Fuse (2 A)		2429
-F3	Fuse (2 A)		2429
	Total of elements 31		

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AZBIL TELSTAR TECHNOLOGIES, S.L.U. guarantees the repair of the defects and the substitution of all the defective pieces IN THE PERIOD OF 2 YEARS, starting at the date of the expedition of the equipment.

This equipment has undergone satisfactorily all the tests and revisions. Therefore, any damage caused by transport, wrong installation or inappropriate use is excluded of this guarantee.

This guarantee does not include any of the consumable materials comprised in the equipment, such as filters, fluids, filaments, etc.

The inappropriate use of the equipment excludes AZBIL TELSTAR TECHNOLOGIES, S.L.U. of any responsibility for eventual direct or indirect damages of any nature.

Likewise, this guarantee will be void if pieces other than the ones supplied by AZBIL TELSTAR TECHNOLOGIES, S.L.U. are used for the repairing of the equipment.


NOTE. This guarantee is valid for all the laboratory freeze-dryers unless otherwise stated. In order to obtain any information about this guarantee, please contact Telstar Technologies Customer Service.

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LYOQUEST ACCESSORIES CONFIGURATION

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Model _____
 Condenser temperature _____
 Voltage _____
 Frequency _____
 Lyologger _____
 Ocean Packing _____
 Vacuum pump _____

X = Void (without this option)

Chamber	Plates: SP3 - 3 Standard HP3 - 3 Heated	Stoppering device: S	Integrated Manifold: M8 - 8 inlets	Larger spacer: LS (optional)	Drip tray: DT (optional)	Chamber order code							
						C220	SP3	X	X	X	DT	(X)	
C220	SP3	-	-	-	-	C220	SP3	X	X	X	DT	(X)	
	HP3 *	-	M8	LS	DT	C220	HP3	X	M8	LS			(X)
	S	-	M8	-	DT	C220	HP3	S	X	-			
		S	M8	-	DT	C220	HP3	S	M8	-			
C220													

Plates type _____
 Stoppering device _____
 Integrated manifold _____
 Larger spacer _____
 Drip tray _____

* Heating function disabled in NO PLUS versions

X = Void (without this option)

Chamber	Plates: HP3 - 3 Heated HP4 - 4 Heated	Stoppering device: S	Larger spacer: LS (optional)	Adapter: A (Required)	Chamber order code						
					C350	HP3	X	LS	A	(X)	
C350	HP3 *	-	LS	A	C350	HP3	X	LS	A	(X)	
	HP4 *	S	-		C350	HP3	S	(X)			
C350	HP3 *	-	-	A	C350	HP4	X	X	(X)		
	HP4 *	S	-		C350	HP4	S	X			
C350											

Plates type _____
 Stoppering device _____
 Larger spacer _____
 Adapter _____

* Heating function disabled in NO PLUS versions

X = Void (without this option)

QUALIFICATION OPTIONS

IQ/OQ
FAT
FAT+ IQ/OQ
SAT
DOCUMENTATION IQ/OQ
WITHOUT QUALIFICATION

WITHOUT QUALIFICATION

M	Adapter: A (Required)	Drip tray: DT (optional)	Manifold order code					
			M8	2xM8	M12	M40		
M8	A	DT	M8	A	DT	(X)		
2xM8			(X)					
M12							(X)	
M40								(X)
2xM40								
M								

Manifold type _____
 Adapter _____
 Drip tray _____

X = Void (without this option)

OTHER ACCESSORIES

Code	Description